

# Cervical Fusion Using PliaFX<sup>®</sup> Prime Moldable Demineralized Fibers

Case performed by: Byron Branch, MD, FAANS

CASE STUDY

Degeneration of the discs and instability in the cervical spine can result in neck pain and radiculopathy causing pain, numbness, or weakness in the chest or arms.<sup>1</sup> If conservative treatment options fail to relieve pain and disability, cervical fusion is a surgical option to address these issues. Autograft harvested from a second surgery site on the patient may be used to support fusion, but it is associated with drawbacks such as donor-site pain, increased infection risk, and increased operating room time.<sup>2</sup> Demineralized allograft bone is an alternative which addresses the risks associated with autograft while still providing an osteoconductive scaffold with osteoinductive potential to support new bone growth and fusion.

LifeNet Health developed PliaFX Prime to provide optimized handling capabilities, undiluted osteoinductive potential, and a hospitable scaffold and void filler for cellular attachment<sup>3-10</sup> — all using 100% allograft bone fibers without synthetic or xenograft carriers. These long, interconnected cortical bone fibers are optimally demineralized using LifeNet Health's patented PAD<sup>®</sup> technology and provide a surface that is rough enough to promote cellular attachment, yet contiguous enough to promote cell spreading and intercellular connection. The superior handling of PliaFX Prime comes from the length and width of the fibers, which are designed to promote malleability, while microhooks allow the fibers to interlock, thereby maintaining the graft's shape and ensuring retention at the implant site.

The following case describes the use of PliaFX Prime in cervical fusion.

## Patient

- 51-year-old female patient with a BMI of 29 kg/m<sup>2</sup>
- Cervical spondylolisthesis, foraminal stenosis, and adult spinal deformity (Figures 1, 2)
- Patient experienced degenerative symptoms for 7 years
- Patient was taking anti-seizure/mood-altering, blood pressure, diuretic, and pain medications concomitantly
- Previous surgical history of ACDF at C5-6

## Procedure

- Two levels (C2 and C3) were treated using the Synthes Symphony<sup>™</sup> Spine System
- 10 cc PliaFX Prime was mixed with approximately 1.1 cc local autograft (9:1)

## Results

- Fusion was assessed radiologically at 6, 12 (Figure 3), and 24 (Figure 4) weeks postoperative

## Conclusion

- All levels were successfully fused at 6 months postoperative (Figure 4)
- There were no notable postoperative complications

# Cervical Fusion Using PliaFX<sup>®</sup> Prime Moldable Demineralized Fibers

CASE STUDY




**Figure 1.**

Preoperative lateral X-ray showing area of interest (red circle) and existing hardware from previous surgery (red arrow).



**Figure 2.**

Preoperative lateral computed tomography (CT) scan showing area of interest (red circle).



**Figure 3.**

Lateral X-ray at 3 months postoperative showing formation of bridging bony callus at the graft implant site (red circle).



**Figure 4.**

Lateral X-ray at 6 months postoperative showing bony fusion at the graft implant site (red circle).

# Cervical Fusion Using PliaFX® Prime Moldable Demineralized Fibers

CASE STUDY

## References

1. Jiang SD, Jiang LS, Dai LY. Degenerative cervical spondylolisthesis: a systematic review. *Int Orthop*. 2011 Jun;35(6):869-75.
2. Heneghan HM, McCabe JP. Use of autologous bone graft in anterior cervical decompression: morbidity & quality of life analysis. *BMC Musculoskelet Disord*. 2009;10:158.
3. Data on file LifeNet Health. ES-17-111-02.
4. Data on file LifeNet Health. ES-17-110.
5. Zhang M, Powers RM, Jr, Wolfinbarger L, Jr. Effect(s) of the demineralization process on the osteoinductivity of demineralized bone matrix. *J Periodontol*. 1997;68(11):1085-1092.
6. Pietrzak WS, Woodell-May J, McDonald N. Assay of bone morphogenetic protein-2, -4, and -7 in human demineralized bone matrix. *J Craniofac Surg*. 2006;17(1):84-90.
7. Turonis JW, McPherson JC, 3rd, Cuenin MF, Hokett SD, Peacock ME, Sharawy M. The effect of residual calcium in decalcified freeze-dried bone allograft in a critical-sized defect in the *Rattus norvegicus* calvarium. *J Oral Implantol*. 2006;32(2):55-62.
8. Herold RW, Pashley DH, Cuenin MF, et al. The effects of varying degrees of allograft decalcification on cultured porcine osteoclast cells. *J Periodontol*. 2002;73(2):213-219.
9. Mott DA, Mailhot J, Cuenin MF, Sharawy M, Borke J. Enhancement of osteoblast proliferation in vitro by selective enrichment of demineralized freeze-dried bone allograft with specific growth factors. *J Oral Implantol*. 2002;28(2):57-66.
10. Murphy MB, Suzuki RK, Sand TT, Chaput CD, Gregory CA. Short Term Culture of Human Mesenchymal Stem Cells with Commercial Osteoconductive Carriers Provides Unique Insights into Biocompatibility. *J Clin Med*. 2013;2(3):49-66.

Results from case studies are not predicative of results in other cases. Results in other cases may vary.  
Please refer to the instructions for use for a complete list of indications, contraindications, warnings, and precautions.

LifeNet Health helps to save lives, restore health and give hope to thousands of patients each year. We are the world's most trusted provider of transplant solutions, from organ procurement to new innovations in bio-implant technologies and cellular therapies—a leader in the field of regenerative medicine, while always honoring the donors and healthcare professionals that allow the healing process.

## LifeNetHealth.org

LifeNet Health, the LifeNet Health logo, PAD, and PliaFX are registered trademarks of LifeNet Health.

©2022 LifeNet Health, Virginia Beach, VA. All rights reserved.

The third party trademarks used herein are the trademarks of their respective owners.  
68-20-370.00